



LAB & DIAGNOSTIC PRODUCTS



Equipment and Testing Supplies for Enhanced Detection of Disease

SIMPLE SOLUTIONS. SIGNIFICANT RESULTS.

We take pride in providing reliable, accurate and durable products that aid in earlier diagnosis and prevention of disease. Review the catalog to see the latest additions to our popular line of diagnostic instruments and equipment, as well as our testing and laboratory products.

Diagnostic Instruments	3-5
Patient Monitors	6-7
Blood Pressure Cuffs & Connectors	8-11
Blood Pressure Monitors & Accessories	12-16
Thermometers	17-19
Stethoscopes	20-22
Defibrillators	23
Electrodes	24-28
Point-of-Care Testing	29-36
Lab Supplies	37-39



Diagnostic Center

Choose the Medline diagnostic center with the components you want.

- » Automatic safety shutdown after 2 minutes
- » All handles come with rheotronic technology
- » Activate by briefly tapping the rotating button to +, immediate 100% light output
- » Continuous variable adjustment by briefly tapping to + or -
- » Diode signals whether handle is on or off

**ASK ABOUT
OUR 30-DAY
RISK-FREE
TRIAL.**

Medline integrated diagnostic system and transformer sets Components chart

Item No.	Wall Transformer	Otosopes		Ophthalmoscopes		Thermometer	Specula Dispenser	Aneroid
	MDS3653300	MDS10565301 L2 LED	MDS10580301 L2 XENON	MDS10571301 L2 XENON	MDS10573301 L3 XENON	MDS3656301	MDS3654	MDSL3655109
MDS36533001	•	•		•				
MDS36533002	•		•	•				
MDS36533003	•	•		•			•	
MDS36533004	•		•	•			•	
MDS36533005	•	•		•			•	•
MDS36533006	•		•	•			•	•
MDS36533007	•	•			•		•	•
MDS36533008	•	•		•		•	•	•
MDS36533009	•	•		•		•	•	•
MDS36533010	•	•			•	•	•	•



MDS3722



MDS3727



MDS2131203



MDS2130203

Diagnostic Sets

Item No.	Description	Pkg.
----------	-------------	------

L2 Otoscope and Ophthalmoscope Diagnostic Sets

MDS3746004	3.5V LED/XL with C handle for Li-batteries/Xenon Opth.	1/ea
MDS3750003	3.5V LED/XL with plug-in handle with ri-accu L/Xenon Opth.	1/ea

L3 Otoscope and Ophthalmoscope Diagnostic Sets

MDS3746006	3.5V LED/XL with C handle for Li-batteries/Xenon Opth.	1/ea
MDS3750004	3.5V LED/XL with plug-in handle with ri-accu L/Xenon Opth.	1/ea

e-scope® Otoscope/Ophthalmoscopes sets (Pocket Instruments)

e-scope ophthalmoscopes are available with vacuum, halogen or LED light sources.

- » Innovative 3.7V LED illumination provides a higher contrast for a more efficient diagnosis
- » Diopter disc with 18 corrective lenses (+/- 1, 2, 3, 4, 6, 8, 10, 15, 20 diopters)
- » Easy-to-use aperture wheel with six different apertures (fixation star, large circle, small circle, red free filter, blue filter, semi circle)
- » High performance optics with aspherical condenser lens
- » Spectacle protection
- » Parallel optical path

Item No.	Description	Pkg.
----------	-------------	------

Fiber optic-2.5V Xenon, Otoscope/Ophthalmoscope

MDS2130202	White, hard case	1/ea
MDS2131202	Black, hard case	1/ea

Fiber optic-3.7V LED, Otoscope/Ophthalmoscope

MDS2130203	White, hard case	1/ea
MDS2131203	Black, hard case	1/ea

Risk-free Trial

Try any Medline diagnostic instrument, set or wall system. We're sure you'll like it, but we'll refund your purchase and pay return shipping if you don't. No restocking fees, no hard feelings. Just return the product in its original packaging.

ri-scope L2 Ophthalmoscope

- » Diopter wheel with 29 corrective lenses:
Plus 1-10, 12, 15, 20, 40, Minus 1-10, 15, 20, 25, 30, 35
- » Easy-to-operate aperture hand-wheel with semi-circle, small/medium/large circle fixation star, slit and grid
- » Optimized high-performance optics with aspherical condenser lens
- » Parallel beam path
- » Dust-proof

Item No.	Description	Pkg.
MDS3725	3.5V XL for Lithium Batteries, C-handle	1/ea

ri-scope L3 Ophthalmoscope

- » Corrective lenses for 89 diopter values:
Plus 1-45 in single steps, Minus 1-44 in single steps
- » Easy-to-operate aperture hand-wheel with semi-circle, small/medium/large circle, fixation star, slit and grid
- » Includes filter wheel engagable for all apertures with symbol display, red-free filter, blue filter and polarization filter
- » Focusing device with zero retention force for the quick correction of ametropias
- » Dust-proof

Item No.	Description	Pkg.
MDS3726	3.5V XL for Lithium Batteries, C-handle	1/ea

ri-scope L2 Otoscope

- » Fiber optics for optimal light transmission
- » Bi-directional swiveling, sealing optical glass with 3x magnification
- » Optional operation lens with 4x magnification available
- » Suitable for pneumatic tests (supplied without ball)

Item No.	Description	Pkg.
MDS3716	2.5V XL, Xenon, for Alkaline Batteries, C-handle	1/ea
MDS3703	3.5V LED for Lithium Batteries, C-handle	1/ea

ri-scope L3 Otoscope

- » Fiber optics for optimal beam and light transmission
- » Bi-directional swiveling optical glass with 3x magnification
- » Operation lens with 4x magnification and glass mini-optic ease the insertion of instruments
- » Easy-to-operate SpecEjec specula ejection device
- » Enlarged conduit for instruments
- » Suitable for pneumatic tests (supplied without ball)

Item No.	Description	Pkg.
MDS3701	2.5V, Xenon, for 2 Alkaline Batteries, C-handle	1/ea
MDS3704	3.5V, LED, for 2 Lithium Batteries, C-handle	1/ea

Speculas

- » Helps reduce cross-contamination

Item No.	Description	Pkg.
MDS14065531	2.5 mm, Disposable, Sleeve	1,000/bg
MDSK10321	2.5 mm	8,500/cs
MDSK10321Z	2.5 mm	850/bx
MDSK10321H	2.5 mm	34/bx
MDS14065534	4.0 mm, Disposable, Sleeve	1,000/bg
MDSK10331	4.25 mm	8,500/cs
MDSK10331Z	4.25 mm	850/bx
MDSK10331H	4.25 mm	34/pk



MDSK10321



MDSK10331

Vital Signs Monitor

Medline's Vital Signs Monitor gives you the ability to provide spot and continuous measurement for NIBP, SpO₂ and temperature in one monitor.

- » Portable and lightweight
- » Easy-to-read and user-friendly monitor
- » Available in different configurations to meet your clinical needs

Choose Your Configuration

Item Number	BP	SpO ₂	Oral Temp	Ear Temp	Printer	Stand	Pkg.
MDSM3ASN1	•	•					1/ea
MDSM3ASNTC	•	•	•				1/ea
MDSM3ASNTCP	•	•	•		•		1/ea
MDSM3ASNTTP	•	•		•	•		1/ea
MDSM3ASN	•	•				•	1/ea
MDSM3ASNT	•	•	•			•	1/ea
MDSM3ASNTP	•	•	•		•	•	1/ea
MDSM3ASNTPIR	•	•		•	•	•	1/ea
MDSM3ASNTT	•	•		•			1/ea
MDSM3ASNTTR	•	•		•		•	1/ea

Choose Your Accessories

Item Number	Description	Pkg.
MDSM3ASTAND01	Rolling Stand w/Storage Basket	1/ea
MDSM3AWALL01	Wall Mount w/Storage Basket	1/ea
MDSSPO2ADL	Adult Pulse Oximeter Finger Sensor, Reusable	1/ea
MDSSPO2DISP	Adult Pulse Oximeter Finger Sensor, Disposable	1/ea
MDSSPO2PED	Pediatric Pulse Oximeter Finger Sensor, Reusable	1/ea
MDSSPO2PEDDSP	Pediatric Pulse Oximeter Finger Sensor, Disposable	1/ea
MDSM3ACUFFI	Edan Blood Pressure Cuff, Infant (14–21 cm)	1/ea
MDSM3ACUFFC	Edan Blood Pressure Cuff, Child (20.5–28 cm)	1/ea
MDSM3ACUFFA	Edan Blood Pressure Cuff, Adult (25–35 cm)	1/ea
MDSM3ACUFFL	Edan Blood Pressure Cuff, Large Adult (33–47 cm)	1/ea
MDSM3ACUFFT	Edan Blood Pressure Cuff, Thigh (42–54 cm)	1/ea
MDSM3AWRRNTY	Additional 1-Year Extended Warranty	1/ea
KDL500500	Filac FasTemp Thermometer Probe Covers	500/cs
KDL500500Z	Filac FasTemp Thermometer Probe Covers	20/bx
MDSM3A7208058	Ear Probe Covers	200/bx
MDS9701	Ear Probe Covers	100/bx

Additional Parts

Item Number	Description	Pkg.
MDSBPTUBE	NIBP Extension Tube, 3 meters	1/ea
MDSSPO2EXT	SpO ₂ Extension Cable, 2 meters	1/ea
MDSM3A101927	Replacement Single Caster Wheel for Vital Signs Stand	1/ea
MDSPOWERCORD	Standard Power Cord for Vital Signs Monitor	1/ea
MDSBATTERY67	Rechargeable Lithium Battery, 14.8V, 4.4AH	1/ea
MDSBATTERY68	Rechargeable Lithium Battery, 14.8V, 2.4AH	1/ea
MDSM3ABASKET	Secondary Basket for Vital Signs Stand	1/ea
MDSM3A78035	Printer Paper	1/ea



MDSM3ASNTP

**SEE PAGES
24-28 FOR OUR
ELECTRODE
OFFERING.**



Medline 12 Channel Interpretive ECG Machine

This 12 channel interpretive unit streamlines patient monitoring with efficient and accurate recordings that will help properly diagnose and monitor patients.

- » Easy-to-read, high resolution color display screen
- » Simplifies patient information entry with one-touch operation and shortcut keys
- » Resists reading contamination by outside noise due to comprehensive filters and high CMRR
- » 2-year warranty
- » Unit moves easily between rooms with the built-in handle or on the optional mobile stand, available separately

Item No.	Description	Pkg.
MDSECGSE1201	Monitor, ECG, 12 Lead, SE 1201	1/ea
MDSECGSTAND01	Stand, ECG, F/SE 1201	1/ea
GCO32017349	Z-Fold ECG Replacement Paper	20/cs
MDS0157107048	Replacement Leadwire Cable	1/ea
MDS0157040173	Replacement Alligator Clips for Leadwire Cable	1/pk
MDSM616101	ECG Resting Tab Electrodes	500/bx



Medline Handheld Fetal and Vascular Doppler

This compact Doppler unit provides important, noninvasive testing to help clinicians diagnose many serious conditions, such as blood clots, heart valve defects, congenital heart disease, blocked or bulging arteries and decreased blood circulation. They are extremely popular tools found in most clinical settings.


- » Noise resistant speakers
- » Lightweight, portable design
- » LCD screen displays FHR and probe frequency
- » Interchangeable waterproof probe design
- » Headphone plug

Item No.	Description	Pkg.
MDSDOPPLERF	Doppler, Fetal, 2 MHZ	1/ea
MDSDOPPLERVF	Doppler, Fetal, 2 MHZ and Vascular 4 MHZ	1/ea
MDSDOPPLERV	Doppler, Vascular, 4 MHZ	1/ea
MDSDOPPLER8	Doppler, Vascular, 8 MHZ	1/ea
MDSPROBEV	Probe, Vascular, 4 MHZ	1/ea
MDSPROBEF	Probe, Fetal, 2 MHZ	1/ea
MDSPROBEV8	Probe, Vascular, 8 MHZ	1/ea
MDSPROBEF3	Probe, Fetal, 3 MHZ	1/ea

Disposable Blood Pressure Cuffs and Connectors

Choosing the right size cuff for each patient ensures accurate readings and allows less room for error.

The chart below will help you determine the correct cuff size based on the patient's upper arm circumference. Please note that if you choose a cuff that is not the appropriate size, readings may not be accurate.



Cuff Size (Limb Circumference)	DINAMAP CONNECTOR (DM)		BAYONET-STYLE CONNECTOR (HP)		SLIP LUER CONNECTOR (SL)	
	1-Tube	2-Tube	1-Tube	2-Tube	1-Tube	2-Tube
Infant (9–14.8 cm)	MDS9910DM	MDS9920DM	MDS9910HP	MDS9920HP	MDS9910SL	MDS9920SL
Infant Long (9–14.8 cm)	Contact Representative	MDS9920DML	MDS9910HPL	Contact Representative	Contact Representative	Contact Representative
Child (13.8–21.5 cm)	MDS9911DM	MDS9921DM	MDS9911HP	MDS9921HP	MDS9911SL	MDS9921SL
Child Long (13.8–21.5 cm)	Contact Representative	MDS9921DML	MDS9911HPL	Contact Representative	Contact Representative	Contact Representative
Small Adult (20.5–28.5 cm)	MDS9912DM	MDS9922DM	MDS9912HP	MDS9922HP	MDS9912SL	MDS9922SL
Adult (27.5–36.5 cm)	MDS9913DM	MDS9923DM	MDS9913HP	MDS9923HP	MDS9913SL	MDS9923SL
Adult Long (27.5–36.5 cm)	MDS9913DML	MDS9923DML	MDS9913HPL	MDS9923HPL	Contact Representative	MDS9923SLL
Large Adult (35.5–46 cm)	MDS9914DM	MDS9924DM	MDS9914HP	MDS9924HP	MDS9914SL	MDS9924SL
Large Adult Long (35.5–46 cm)	MDS9914DML	MDS9924DML	MDS9914HPL	MDS9924HPL	MDS9914SLL	MDS9924SLL
Thigh (45–56.5 cm)	MDS9915DM	MDS9925DM	MDS9915HP	MDS9925HP	MDS9915SL	MDS9925SL

Cuffs constructed of vinyl are available on items listed above. Place a “V” after the connector type code to order. Yellow vinyl cuffs are available in select sizes for isolation areas. Contact your Medline representative for more information.



MALE LUER CONNECTOR (ML)		MARQUETTE FEMALE CONNECTOR (MQ)		MARQUETTE MALE CONNECTOR (MQML)	MARQUETTE MALE & FEMALE CONNECTOR (MQMF)
Recommended for: CasMed® 740/750 units		Recommended for: GE (Marquette®) and Welch Allyn units		Recommended for: older Spacelabs (Instrumentarium), Datascope, Dräger, Fukuda Denshi, BCI, Colin, Welch Allyn and Nihon Kohden units	Recommended for: Datex-Ohmeda® units
1-Tube	2-Tube	1-Tube	2-Tube	1-Tube	2-Tube
MDS9910ML	MDS9920ML	MDS9910MQ	MDS9920MQ	MDS9910MQML	MDS9920MQMF
Contact Representative	Contact Representative	Contact Representative	MDS9920MQL	Contact Representative	Contact Representative
MDS9911ML	MDS9921ML	MDS9911MQ	MDS9921MQ	MDS9911MQML	MDS9921MQMF
Contact Representative	Contact Representative	Contact Representative	MDS9921MQL	Contact Representative	Contact Representative
MDS9912ML	MDS9922ML	MDS9912MQ	MDS9922MQ	MDS9912MQML	MDS9922MQMF
MDS9913ML	MDS9923ML	MDS9913MQ	MDS9923MQ	MDS9913MQML	MDS9923MQMF
MDS9913MLL	MDS9923MLL	MDS9913MQL	MDS9923MQL	MDS9913MQMLL	MDS9923MQMFL
MDS9914ML	MDS9924ML	MDS9914MQ	MDS9924MQ	MDS9914MQML	MDS9924MQMF
MDS9914MLL	MDS9924MLL	MDS9914MQL	MDS9924MQL	MDS9914MQMLL	MDS9924MQMFL
MDS9915ML	MDS9925ML	MDS9915MQ	MDS9925MQ	MDS9915MQML	MDS9925MQMF

Blood Pressure Cuff Covers and Barriers

Help reduce the potential of cross-contamination by using protective, disposable blood pressure cuff covers.

- » Manufactured without latex
- » Made of Tyvek® material
- » Quick and easy to replace or remove

Item No.	Description	Pkg.
Cuff Covers		
MDS9160	Child Cuff Cover	20/bx
MDS9161	Adult Cuff Cover	20/bx
MDS9162	Large Adult Cuff Cover	20/bx
Cuff Barriers		
MDS917113	Small Cuff Barrier	300/cs
MDS917114	Medium Cuff Barrier	300/cs
MDS917115	Large Cuff Barrier	150/cs









MDS9161

Adapters

Streamline the conversion process with the correct adapter.

- » Metal adapters are more durable than plastic
- » Permanently heat-sealed and replaced free of charge
- » In-service staff training for adapter installation

<p>ITEM NUMBER MDS8230</p> <p>Metal Bayonet Adapter 5/32" barbed single tube to female HP connection</p> 	<p>ITEM NUMBER MDS182YRX4</p> <p>Y Tube Bayonet Adapter 1/8" barbed double tubes to female HP connection</p> 	<p>ITEM NUMBER STCZFTBPF240</p> <p>Plastic Bayonet Adapter 5/32" single tube, MRI safe</p> 
<p>ITEM NUMBER MDS8220</p> <p>Metal Bayonet Adapter 1/8" barbed single tube to female HP connection</p> 	<p>ITEM NUMBER MDSDM2YRX4</p> <p>Y Tube DM to Bayonet Adapter Male DM double tubes to female HP connection</p> 	<p>ITEM NUMBER STCAD182YBPF</p> <p>Plastic Bayonet Adapter 1/8" barbed double tube to female HP connection, MRI safe</p> 

Reusable Blood Pressure Cuffs

Reuse with Confidence

We build our reusable cuffs with a strong polymer that can withstand up to 10,000 inflations. Then we add rounded edges for patient comfort and color coding for easy size identification. The result? An intuitive, comfortable cuff that you can depend on time and time again.

The Right Size Every Time

Choosing the right size cuff for each patient ensures accurate readings and allows less room for error. The chart below will help you determine the correct cuff size based on the patient's upper arm circumference. Please note that if you choose a cuff that is not the appropriate size, readings may not be accurate.



	Bayonet-Style Connector	Screw Connector	
	1-Tube HP	1-Tube DM	2-Tube DM
Infant (9–14.8 cm)	MDS9910HPRE	MDS9910DMRE	MDS9920DMRE
Infant Long (9–14.8 cm)	MDS9910HPLRE	MDS9910DMLRE	MDS9920DMLRE
Child (13.8–21.5 cm)	MDS9911HPRE	MDS9911DMRE	MDS9921DMRE
Child Long (13.8–21.5 cm)	MDS9911HPLRE	MDS9911DMLRE	MDS9921DMLRE
Small Adult (20.5–28.5 cm)	MDS9912HPRE	MDS9912DMRE	MDS9922DMRE
Small Adult Long (20.5–28.5 cm)	MDS9912HPLRE	Contact Representative	MDS9922DMLRE
Adult (27.5–36.5 cm)	MDS9913HPRE	MDS9913DMRE	MDS9923DMRE
Adult Long (27.5–36.5 cm)	MDS9913HPLRE	MDS9913DMLRE	MDS9923DMLRE
Large Adult (35.5–46 cm)	MDS9914HPRE	MDS9914DMRE	MDS9924DMRE
Large Adult Long (35.5–46 cm)	MDS9914HPLRE	MDS9914DMLRE	MDS9924DMLRE
Thigh (45–56.5 cm)	MDS9915HPRE	MDS9915DMRE	MDS9925DMRE

	Slip Luer Connector	Marquette Connector	
	1-Tube SL	1-Tube MQRE	2-Tube MQRE
Infant (9–14.8 cm)	MDS9910SLRE	MDS9910MQRE	MDS9920MQRE
Infant Long (9–14.8 cm)	MDS9910SLLRE	MDS9910MQLRE	MDS9920MQLRE
Child (13.8–21.5 cm)	MDS9911SLRE	MDS9911MQRE	MDS9921MQRE
Child Long (13.8–21.5 cm)	MDS9911SLLRE	MDS9911MQLRE	MDS9921MQLRE
Small Adult (20.5–28.5 cm)	MDS9912SLRE	MDS9912MQRE	MDS9922MQRE
Small Adult Long (20.5–28.5 cm)	Contact Representative	Contact Representative	MDS9922MQLRE
Adult (27.5–36.5 cm)	MDS9913SLRE	MDS9913MQRE	MDS9923MQRE
Adult Long (27.5–36.5 cm)	MDS9913SLLRE	MDS9913MQLRE	MDS9923MQLRE
Large Adult (35.5–46 cm)	MDS9914SLRE	MDS9914MQRE	MDS9924MQRE
Large Adult Long (35.5–46 cm)	MDS9914SLLRE	MDS9914MQLRE	MDS9924MQLRE
Thigh (45–56.5 cm)	MDS9915SLRE	MDS9915MQRE	MDS9925MQRE

medline**elite**

Automatic Digital Blood Pressure Monitor

Revolutionary design and proprietary programming provide automatic inflation and error detection.

- » Customized inflation levels for each individual patient
- » Color gauge for quick detection
- » Provides quick digital readings of systolic/diastolic pressure and pulse rate
- » Easy-to-use one-push button design
- » Stores the 90 most recent measurements
- » Large digital readout
- » Includes the automatic blood pressure unit, cuff and four AA batteries

Item No.	Description	Pkg.
MDS3001	Monitor with Adult Cuff	1/ea
MDS3001LA	Monitor with Large Adult Cuff	1/ea
MDS3001PLUS	Monitor with Adult and Large Adult Cuffs	1/ea
MDS4001	Monitor with Adult Cuff	1/ea
MDS4001LAT	Talking Monitor with Large Adult Cuff	1/ea
MDS9970	Small Adult Cuff (fits arms 17–22 cm)	1/ea
MDS9971	Adult Cuff (fits arms 22–30 cm)	1/ea
MDS9972	Large Adult Cuff (fits arms 30–42 cm)	1/ea
MDS9973	Extra Large Adult Cuff (fits arms 42–48 cm)	1/ea
MDS9975	Universal Cuff (fits arms 22–42 cm)	1/ea

Easy Product Selection

Differences among various aneroids can be subtle. The names of some of our popular aneroid products reflect our good, better and best in each line.

medline**pro**

Semi-Automatic Digital Blood Pressure Monitor

Reliable and accurate readings from a manually inflated digital blood pressure unit.

- » Provides quick digital readings of systolic/diastolic pressure and pulse rate
- » Easy-to-use one-push button design
- » Stores the 90 most recent measurements
- » Unit includes the semi-automatic blood pressure unit, cuff and four AA batteries

Item No.	Description	Pkg.
MDS3002	Monitor with Adult Cuff	1/ea
MDS3002LA	Monitor with Large Adult Cuff	1/ea
MDS9770	Small Adult Cuff (fits arms 17–22 cm)	1/ea
MDS9771	Adult Cuff (fits arms 22–30 cm)	1/ea
MDS9772	Large Adult Cuff (fits arms 30–42 cm)	1/ea
MDS9773	Extra Large Adult Cuff (fits arms 42–48 cm)	1/ea

medline**plus**

Digital Wrist Blood Pressure Monitor

Easy-to-use automatic wrist blood pressure monitor.

- » Automatically inflates and deflates
- » Provides quick digital readings of systolic/diastolic pressure and pulse rate
- » Memory with date and time
- » Color gauge for quick detection
- » Includes two AAA batteries

Item No.	Description	Pkg.
MDS3003	Digital Wrist Blood Pressure Monitor	1/ea



MDS3001PLUS



MDS4001



MDS3003



MDS9410

medline**elite** Medline Handheld Aneroids

Our finest handheld aneroids.

- » Luminescent dial is easy to read, even in dim light
- » Shatter-resistant acrylic faceplate surrounds the luminescent dial
- » Exceptional accuracy from a top-quality, zinc alloy precision gauge
- » Nylon range finder cuff for quick, convenient patient fitting
- » Carrying case included

Item No.	Description	Pkg.
MDS9411	Medline Elite Handheld Aneroid, Child	1/ea
MDS9410	Medline Elite Handheld Aneroid, Adult	1/ea
MDS9413	Medline Elite Handheld Aneroid, Large Adult	1/ea



MDS9380LF

medline**pro** Neoprene Handheld Aneroids

High-quality handheld aneroids.

- » Durable neoprene tubing and bladder
- » Nylon range finder cuff for quick, convenient patient fitting
- » Carrying case included

Item No.	Description	Pkg.
MDS9386LF	Medline Pro Handheld Aneroid, Infant	1/ea
MDS9387LF	Medline Pro Handheld Aneroid, Child	1/ea
MDS9380LF	Medline Pro Handheld Aneroid, Adult	1/ea
MDS9388LF	Medline Pro Handheld Aneroid, Large Adult	1/ea



MDS9380

medline**plus** PVC Handheld Aneroids

Great quality. Great price.

- » PVC tubing and bladder
- » Nylon range finder cuff for quick, convenient patient fitting
- » Carrying case included

Item No.	Description	Pkg.
MDS9386	Medline Plus Handheld Aneroid, Infant	1/ea
MDS9387	Medline Plus Handheld Aneroid, Child	1/ea
MDS9380	Medline Plus Handheld Aneroid, Adult	1/ea
MDS9388	Medline Plus Handheld Aneroid, Large Adult	1/ea

Wall-Mount/Mobile Aneroids

Large, raised numbers and a luminescent coating make the wall-mount/mobile aneroid easy to read in dim light.

- » Eight-foot coiled tubing and pivoting bracket extend service area
- » Inflation system and bulb assembly store easily and safely behind gauge
- » PVC and Neoprene are manufactured without latex
- » Includes Medline's nylon range finder cuff

Item No.	Description	Pkg.
medlinepro Neoprene Aneroids		
MDS9400LF	Aneroid, Wall Mount, Adult Cuff	1/ea
MDS9407LF	Aneroid, Mobile, Adult Cuff with Stand	1/ea
MDS9418	Mobile Stand	1/ea
medlineplus PVC Aneroids		
MDS9400	Aneroid, Wall Mount, Adult Cuff	1/ea
MDS9407	Aneroid, Mobile, Adult Cuff with Stand	1/ea
MDS9418	Mobile Stand	1/ea



MDS9400LF

medlineplus

Home Blood Pressure Kits

- » Easy-to-read aneroid gauge
- » Includes complete instructions and a recording chart
- » Packed in a matching nylon carrying case

Item No.	Description	Pkg.
MDS9300	Home Blood Pressure Kit, Detached Stethoscope	1/ea
MDS9301	Home Blood Pressure Kit, Attached Stethoscope	1/ea

medlineplus

Palm Aneroids

Hand-held palm aneroid swivels for both left and right hand users

Item No.	Description	Pkg.
MDSANERIOD1	Handheld Palm Aneroid	1/ea
MDSANERIOD10	Child Cuff	1/ea
MDSANERIOD11	Adult Cuff	1/ea
MDSANERIOD12	Large Adult Cuff	1/ea



MDSANERIOD1

medline**elite**

Compli-Mates® Combination Kits

Kit includes an aneroid sphygmomanometer, matching 22" (56 cm) Sprague-Rappaport stethoscope and matching nylon carrying case.

Item No.	Description	Pkg.
MDS9121	Magenta	1/ea
MDS9123	Purple	1/ea
MDS9124	Royal Blue	1/ea
MDS9125	Black	1/ea
MDS9128	Hunter Green	1/ea
MDS9140	Steel/Black (Stainless Steel)	1/ea



MDS9124

medline**pro**

Compli-Mates Combination Kits

Kit includes an aneroid sphygmomanometer, matching dual-head stethoscope and matching nylon carrying case.

Item No.	Description	Pkg.
MDS9116	Magenta	1/ea
MDS9117	Teal	1/ea
MDS9118	Purple	1/ea
MDS9119	Royal Blue	1/ea
MDS9120	Black	1/ea
MDS9133	Hunter Green	1/ea
MDS9136	Red	1/ea



MDS9117

medline**plus**

Compli-Mates Aneroid Sphygmomanometers

Kit includes an aneroid sphygmomanometer and matching nylon carrying case.

Item No.	Description	Pkg.
MDS9104	Red	1/ea
MDS9106	Navy Blue	1/ea
MDS9108	Hunter Green	1/ea
MDS9109	Orchid	1/ea
MDS9111	Magenta	1/ea
MDS9112	Teal	1/ea
MDS9113	Purple	1/ea
MDS9114	Royal Blue	1/ea



MDS9104

Aneroid Accessories



MDS9414



MDS9630



MDS9508

Handheld Aneroids			Wall-Mount Aneroids	
medline elite	medline pro	medline plus	medline pro	medline plus
MDS9410, MDS9411, MDS9413	MDS9386LF, MDS9387LF, MDS9380LF, MDS9388LF	MDS9386, MDS9387, MDS9380, MDS9388	MDS9407LF, MDS9400LF	MDS9400, MDS9407

Premier Gauge and Parts

Gauge (MDS9414)	•				
Glass Face (MDS9415)	•				

Standard Gauge and Parts

Aneroid Gauge (MDS9381)		•	•		
Glass Face (MDS9507)		•	•		

Bulbs and Valves

Valve (MDS9508)	•	•	•	•	•
PVC Bulb (MDS91481)			•		•
Neoprene Bulb (MDS91481LF)		•		•	
Premier bulb and valve assembly (MDS9630)	•				
PVC Bulb and Valve Assembly (MDS91482)		•	•	•	•

Two-Tube Inflation Bag and Nylon Range Finder Cuff

Neoprene, Infant (MDS91422LF)		•		•	
Neoprene, Child (MDS91421LF)		•		•	
Neoprene, Adult (MDS91420LF)		•		•	
Neoprene, Large Adult (MDS91423LF)		•		•	
PVC, Infant (MDS91422)			•		•
PVC, Child (MDS91421)			•		•
PVC, Adult (MDS91420)			•		•
PVC, Large Adult (MDS91423)			•		•
PVC, Thigh (MDS91424)			•		•

Two-Tube Seamless Inflation Bag

Neoprene, Adult (MDS9505LF)		•		•	
Neoprene, Large Adult (MDS9506LF)		•		•	
PVC, Adult (MDS9505)			•		•
PVC, Large Adult (MDS9506)			•		•

Nylon Range Finder Hook-and-Loop Cuffs

Child (MDS9421)		•	•	•	•
Adult (MDS9420)		•	•	•	•
Large Adult (MDS9423)		•	•	•	•

8' Coiled Tubing with Connector

PVC (MDS9475)					•
Neoprene (MDS9475LF)				•	

Miscellaneous Parts

Metal Connector/Coupler (MDS9481)				•	•
Basket for Wall-Mounted or Mobile Aneroids (MDS9472)				•	•
Mobile Stand				•	•

No Touch Forehead Thermometer

Innovative, quick-response ear thermometer features:

- » One-touch operation
- » Sensor motion readings at 1-2 cm when parallel to skin
- » Reduces cross contamination by eliminating patient contact
- » No probe covers needed
- » Displays green light for healthy temperature and red light for fever
- » Has human mode (child symbol) and surface mode (thermometer symbol)

Item No.	Description	Pkg.
MDSTEMPSCAN	Forehead Thermometer	1/ea



MDSTEMPSCAN

Digital Temple Thermometer

The Medline instant read digital temple thermometer accurately measures body heat using R.A.T.E.™ technology (Rapid, Accurate, Temperature Establishment).

- » Sensors rapidly track heat flow generated from blood vessels to the skin's surface and convert measurements to body temperature
- » Easy-to-use on asleep or alert patients
- » Meets accuracy standards
- » Large digital display
- » Uses two AAA batteries
- » Memory stores previous reading

Item No.	Description	Pkg.
MDS9698L	Temple Thermometer	1/ea
MDS9699	Probe Covers	250/bx



MDS9698



MDS9700

Tympanic Thermometer

Innovative, quick-response ear thermometer with quick probe release.

- » One-touch operation
- » Readings in six to ten seconds depending on individual
- » Mechanical probe cover quick-release allows easy probe cover change
- » Safe to use on patients of any age
- » Readings in Fahrenheit and Celsius
- » Takes an accurate, non-invasive temperature from the tympanic membrane

Item No.	Description	Pkg.
MDS9700	Thermometer	1/ea
MDS9701	Probe covers	100/bx



EXGTAT2000



EXGTAT5000

Exergen® Temporal Scanner™ Thermometer

Precision infrared scanner provides a truly non-invasive method for measuring temperature. Simply scan across the temporal artery area of the forehead, touch behind the earlobe, and an accurate reproducible temperature is measured in about three seconds.

- » Convenient one-handed operation
- » Completely non-invasive
- » Stainless steel probe for highest durability
- » Easy snap-on probe caps
- » Three-year warranty for EXGTAT2000

Item No.	Description	Pkg.
EXGTAT5000	Infrared Temporal Artery Thermometer	1/ea
EXGTAT2000	Infrared Temporal Artery Thermometer	1/ea
EXG134203	Disposable Probe Caps	1,000/cs
EXG134203H	Disposable Probe Caps for EXGTAT2000/5000	50/pk

30-Second Digital Thermometers

Our digital thermometers provide accurate temperatures in about 30 seconds, depending on model.

- » Large LCD screen for easy readings
- » Simple push-button start
- » Disposable sheaths allow for sanitary oral, rectal or axillary use

Item No.	Description	Pkg.
Oral		
MDS9950	Oral °F & °C, Plastic Cover	144/cs
MDS9950Z	Oral °F & °C, Plastic Cover	12/bx
MDS9950H	Oral °F & °C, Plastic Cover	1/ea
MDS9950B	Oral °F & °C, Bulk	200/cs
MDS9950BZ	Oral °F & °C, Bulk	20/bx
MDS9607	Oral Sheaths	20/bx
MDS9928	Oral Kit w/20 Sheaths, °F & °C	1/ea

Rectal

MDS9952	Rectal °F & °C, Plastic Cover	144/cs
MDS9952Z	Rectal °F & °C, Plastic Cover	12/bx
MDS9952H	Rectal °F & °C, Plastic Cover	1/ea
MDS9952B	Rectal °F & °C Bulk	200/cs
MDS9952BZ	Rectal °F & °C Bulk	20/bx
MDS9608	Rectal Sheaths with Lube	100/bx
MDS9929	Rectal Kit w/20 Sheaths, °F & °C	1/ea

Flex-tip

MDS99901	Flex-tip, Waterproof, °F & °C	240/cs
MDS99901Z	Flex-tip, Waterproof °F & °C	12/bx
MDS99901H	Flex-tip, Waterproof, °F & °C	1/ea



MDS9950



MDS9952



MDS99901

Disposable Penlight

Lightweight and disposable, these high-intensity exam lights are perfect for the physician, nurse and EMT. They're constructed of high-impact molded styrene and can be clipped on pockets.

Item No.	Description	Pkg.
MDS131040	Disposable Penlight, White	6/pk
MDS131040H	Disposable Penlight, White	1/ea



MDS131040



MDS926501

Synergy Classic® Dual Frequency Stethoscope

Our Synergy Classic Dual Frequency Stethoscope offers an elevated diaphragm that provides excellent auscultation and acoustics. Apply minimal pressure to hear low frequencies and increase pressure for higher frequencies.

- » High end stainless steel binaurals and chest piece help trap sound
- » Durable, thick wall tubing to eliminate sound interference
- » Soft, narrow ear tips for everyday comfort

Item No.	Description	Pkg.
MDS926501	Synergy Classic Dual Frequency Stethoscope	1/ea



MDS926203

Medline Plus Value Stethoscopes

Quality and economy make our single and dual-head stethoscopes an excellent value.

- » Both stethoscopes feature 22" (56 cm) tubing and lightweight design
- » Thick-walled tubing minimizes background interference
- » Dual-heads feature a non-chill ring on the bell for increased patient comfort
- » Manufactured without latex

Item No.	Description	Pkg.
----------	-------------	------

Single-head

MDS926101	Black	1/ea
MDS926102	Gray	1/ea
MDS926103	Blue	1/ea
MDS926105	Lavender	1/ea
MDS926106	Red	1/ea
MDS926107	Hunter Green	1/ea
MDS926108	Yellow	1/ea
MDS926109	Pink	1/ea

Dual-head

MDS926201	Black	1/ea
MDS926202	Gray	1/ea
MDS926203	Blue	1/ea
MDS926205	Lavender	1/ea
MDS926206	Red	1/ea
MDS926207	Hunter Green	1/ea
MDS926208	Yellow	1/ea
MDS926209	Pink	1/ea

Medline Pro Stethoscope

The Sprague-Rappaport stethoscope incorporates versatility into an excellent sound-transmitting device.

- » Each stethoscope comes complete with a component kit that includes two sets of ear tips, large and small diaphragm, adult bell, pediatric bell, and infant bell
- » The dual 22" (56 cm) tubing design transmits even the faintest sounds

Item No.	Description	Pkg.
MDS926301	Black	1/ea
MDS926302	Gray	1/ea
MDS926303	Blue	1/ea
MDS926305	Lavender	1/ea
MDS926306	Red	1/ea
MDS926307	Hunter	1/ea
MDS926311	Yellow	1/ea
MDS9552	Single/Dual Head Replacement Diaphragm	1/ea



MDS926303

Medline Elite Stainless Steel Stethoscopes

Each stethoscope is handcrafted from high tensile strength stainless steel. Features durable 22" (56 cm) tubing.

- » Color coordinated non-chill rings provide excellent patient comfort
- » All models come in a foam-lined box with one extra set of ear tips
- » Available in adult, pediatric, neonatal and cardiology styles
- » Cardiology stethoscope has a deep, cone-shaped bell – a unique feature that provides excellent sound transmission and acoustic sensitivity. Features 17" (44 cm) tubing

Item No.	Description	Pkg.
Adult		
MDS92250	Gray	1/ea
MDS92260	Black	1/ea
MDS92270	Blue	1/ea
MDS92280	Burgundy	1/ea
MDS92290	Pink Breast Cancer Awareness	1/ea
Pediatric		
MDS92200	Gray	1/ea
Neonatal		
MDS92150	Gray	1/ea
MDS92160	Pink	1/ea
Cardiology		
MDS92500	Black	1/ea



MDS92500





MDS9515

Bowles Stethoscope

Our Bowles Stethoscope features thick-walled 22" (56 cm) tubing, a chrome-plated binaural and a chrome-plated, raised-metal-stem chest piece.

Item No.	Description	Pkg.
MDS9515	Black, Adult	1/ea

Pediatric/Neonatal Stethoscope

The pediatric and neonatal stethoscopes have all the features of the dual-head stethoscopes, but are sized to fit the needs of the individual practitioner and patient.

- » Both stethoscopes feature 22" (56 cm) tubing and lightweight design
- » Thick-walled tubing minimizes background interference
- » Aluminum binaurals

Item No.	Description	Pkg.
MDS9557	Gray, Pediatric	1/ea
MDS9565	Gray, Neonatal	1/ea




MDS9557

Disposable Stethoscope

- » Designed for single patient use
- » Helps reduce cross-contamination in infectious areas
- » High visibility yellow

Item No.	Description	Pkg.
MDS9543	Disposable Stethoscope	1/ea



MDS9543



DFTDCF1310EN



DFTDCFA120



DFTDCFA100EN

Defibtech® Lifeline Series

Built to the demanding standards of professional emergency personnel, Lifeline AEDs also meet US military “drop and shock” specifications.

Item No.	Description	Pkg.
DFTDCF1310EN	Lifeline VIEW® AED (Semi-Automatic)	1/ea
DFTDCFA100EN	Lifeline AED® (Semi-Automatic)	1/ea
DFTDCFA120	Lifeline AUTO AED (Fully-Automatic)	1/ea



DFTDDP100



DFTDD200P

Defibtech Replacement Accessories

Item No.	Description
DFTDDP100	Lifeline AED Adult Pads for adult use only (eight years or older)
DFTDDP200P	Lifeline AED Pediatric Pads for infants and children under eight years old
DFTDDP2001	Lifeline VIEW AED Adult Pads for adult use only (eight years or older)
DFTDDP2002	Lifeline VIEW AED Pediatric Pads for infants and children under eight years old
DFTDCF200	Lifeline AED Battery Pack This replacement battery pack (DBP-1400) has a standby life (installed in the unit) of five years. Each battery pack comes with a new DAC-410 9-volt lithium battery.
DFTDCF2003	Lifeline VIEW AED Battery Pack This replacement battery pack has a standby life (installed in the unit) of four years.



DFTDDP2001



DFTDDP2002



DFTDCF200



DFTDCF2003

General Monitoring Electrodes

- » Clear tracing every time with our high performance, specially formulated gel with extended surface area
- » Electrodes stay securely on the skin, yet remove cleanly with our top-quality adhesive
- » Flexible lead wire attachment with our perforated strip backing
- » Prolonged electrode freshness with our resealable bag (50 pack only)
- » Easy removal with full-edge lift tab

MedGel® General Monitoring Cloth Electrode

Recommended for: ICU/CCU, Emergency, O.R., Patient Rooms
Great for patients with sensitive skin.

- » Enhanced comfort by conforming to the shape of the body with our non-woven breathable cloth material

Item No.	Description	Pkg.
MDSM611303	General Monitoring, Cloth	3/pk, 60/bx, 600/ct, 3,600/cs
MDSM611305	General Monitoring, Cloth	5/pk, 100/bx, 1,000/ct, 6,000/cs
MDSM611350	General Monitoring, Cloth	50/pk, 1,000/ct, 6,000/cs



MDSM611350

MedGel General Monitoring Foam Electrode

Recommended for: ICU/CCU, Emergency, O.R., Patient Rooms
Our best-selling electrode.

- » Flexible, fluid-resistant foam material

Item No.	Description	Pkg.
MDSM611503	General Monitoring, Foam	3/pk, 60/bx, 600/ct, 3,600/cs
MDSM611505	General Monitoring, Foam	5/pk, 100/bx, 1,000/ct, 6,000/cs
MDSM611550	General Monitoring, Foam	50/pk, 1,000/ct, 6,000/cs



MDSM611550

MedGel General Monitoring Tape Electrode

Recommended for: ICU/CCU, Emergency, O.R., Patient Rooms
Ideal for seeing directly under the electrode without having to remove it.

- » Provides visibility under the electrode with our transparent breathable tape material

Item No.	Description	Pkg.
MDSM611403	General Monitoring, Tape	3/pk, 60/bx, 600/ct, 3,600/cs
MDSM611405	General Monitoring, Tape	5/pk, 100/bx, 1,000/ct, 6,000/cs
MDSM611450	General Monitoring, Tape	50/pk, 1,000/ct, 6,000/cs



MDSM611450

Radiotranslucent Electrodes

- » Can be used during X-rays and MRIs in areas where metal is prohibited—Radiology, Cath Lab, Surgery, etc.—with our radiotranslucent snap
- » Clear tracings every time with our high performance, specially formulated gel and extended surface area
- » Stays securely on the skin yet removes cleanly with our top-quality adhesive
- » Prolonged electrode freshness and life with our resealable bag (50 pack only)
- » Flexible lead wire attachment with our perforated strip backing
- » Easy removal with full-edge lift tab

MedGel Radiotranslucent Cloth Electrode

Recommended for: Cath Lab, Radiology
Great for patients with sensitive skin.

- » Enhanced comfort by conforming to the shape of the body with our non-woven breathable cloth material

Item No.	Description	Pkg.
MDSM611250R	Radiotranslucent Cloth	50/pk, 1,000/ct, 6,000/cs



MDSM611250R

MedGel Radiotranslucent Foam Electrode

Recommended for: Cath Lab, Radiology
Our most popular radiotranslucent electrode.

- » Flexible fluid-resistant foam material

Item No.	Description	Pkg.
MDSM611803R	Radiotranslucent Foam	3/pk, 60/bx, 600/ct, 3,600/cs
MDSM611805R	Radiotranslucent Foam	5/pk, 100/bx, 1,000/ct, 6,000/cs
MDSM611850R	Radiotranslucent Foam	50/pk, 1,000/ct, 6,000/cs



MDSM611850R

MedGel Radiotranslucent Tape Electrode

Recommended for: Cath Lab, Radiology
Ideal for seeing directly under the electrode without having to remove it.

- » Provides visibility under the electrode with our transparent breathable tape material

Item No.	Description	Pkg.
MDSM611650R	Radiotranslucent Tape	50/pk, 1,000/ct, 6,000/cs



MDSM611650R

MedGel Pediatric Foam Electrode

Recommended for: PICU

Smaller footprint fits pediatric patients better.

- » Greater percentage of gel allows for gentle use and removal
- » Clear tracings every time with our high performance, specially formulated gel
- » Prolonged electrode freshness and life with our resealable bag
- » Flexible lead wire attachment with our perforated strip backing
- » Easy removal with lift tab

Item No.	Description	Pkg.
MDSM611050	Pediatric, Foam	50/pkg, 1,000/ct, 6,000/cs



MDSM611050


MedGel Pediatric Radiotranslucent Electrode

Recommended for: Cath Lab, Radiology, PICU

Radiotranslucent foam and small size for pediatric patients.

- » Can be used during X-Rays and MRIs in areas where metal is prohibited—Radiology, Cath Lab, Surgery, etc.—with our radiotranslucent snap
- » Greater percentage of gel allows for gentle use and removal
- » Clear tracings every time with our high performance, specially formulated gel
- » Prolonged electrode freshness and life with our resealable bag
- » Flexible lead wire attachment with our perforated strip backing
- » Easy removal with lift tab

Item No.	Description	Pkg.
MDSM6112250R	Pediatric, Radiotranslucent Foam	50/pkg, 1,000/ct, 6,000/cs



MDSM6112250R

MedGel Pre-Wired Neonatal Electrodes

Recommended for: NICU

Pre-wired baby-sized electrodes.

- » Gentle and breathable tape material—ideal for neonatal patients
- » Painless removal with 100% gel coverage
- » Three electrodes come pre-wired for convenience
- » Clear tracings every time with our high performance, specially formulated gel

Item No.	Description	Pkg.
MDSM61123	Neonatal Clear Tape	3/pkg, 300/ct, 3,600/cs
MDSM61143R	Neonatal, Radiotranslucent Clear Tape	3/pkg, 300/ct, 3,600/cs



MDSM61143R

MedGel EKG Resting Tab Electrodes

Recommended for: Cardiology

- » Stable, reliable EKG signals with silver/silver chloride gel
- » Maximum contact with the skin from a large solid conductive area

Item No.	Description	Pkg.
MDSM616101	Tab Electrode, Small	100/pkg, 500/bx, 5,000/cs
MDSM616201	Tab Electrode, Adult	100/pkg, 500/bx, 5,000/cs



MDSM616101




MDSM616201

MedGel Stress/Diaphoretic Foam Electrodes

Recommended for: Cardiology, Stress Test

- » Specially designed for situations where a strong adhesive is required
- » Clear tracings every time with our high performance, specially formulated gel and extended surface area

Item No.	Description	Pkg.
MDSM611903	ECG Stress Foam	3/pk, 60/bx, 600/ct, 3,600/cs
MDSM611905	ECG Stress Foam	5/pk, 100/bx, 1,000/ct, 6,000/cs
MDSM611950	ECG Stress Foam	50/pk, 1,000/bx, 6,000/cs




MDSM611950

MedGel Radiotranslucent Stress/Diaphoretic Foam Electrodes

Recommended for: Cardiology, Cath Lab, Radiology, Stress Test

- » Radiotranslucent and designed for situations where a strong adhesive is required.
- » Can be used during X-Rays and MRIs in areas where metal is prohibited—Radiology, Cath Lab, Surgery, etc.—with our radiotranslucent snap.
- » Clear tracings every time with our high performance, specially formulated gel and extended surface area.

Item No.	Description	Pkg.
MDSM611703R	Radiotranslucent Stress Foam	3/pk, 60/bx, 600/ct, 3,600/cs
MDSM611705R	Radiotranslucent Stress Foam	5/pk, 100/bx, 1,000/ct, 6,000/cs
MDSM611750R	Radiotranslucent Stress Foam	50/pk, 1,000/bx, 6,000/cs



MDSM611750R

MedGel Wet Gel Stress/Diaphoretic Foam Electrodes

Recommended for: Cardiology, Stress Test, EMS

- » Quick reading and stabilization for situations of particularly heavy patient sweating.
- » Instant, high quality trace pickup with our liquid gel

Item No.	Description	Pkg.
MDSM618303	Wet Gel Stress Foam	3/pk, 60/bx, 600/ct, 3,600/cs
MDSM618305	Wet Gel Stress Foam	5/pk, 100/bx, 1,000/ct, 6,000/cs
MDSM618350	Wet Gel Stress Foam	50/pk, 1,000/bx, 6,000/cs



MDSM618350

Medline Plus Repositionable Electrodes

Recommended for: ICU/CCU, Emergency, O.R., Cath Lab, Radiology, PICU, Patient Rooms

- » Easy to reposition at least once with a conductive hydrogel
- » Gentle on sensitive skin
- » Carbon snap for standardization throughout your facility













Item No.	Description	Pkg.
MDS611300	Adult, Solid Gel	5/pk, 100/bx, 1,000/cs



MDS611300

MedGel Electrodes

Recommended Electrodes for Optimal Utilization

		ICU/CCU	Emergency	O.R.	Cardiology	Cath Lab	Radiology	Stress Test	PICU	NICU	Patient Rooms	EMS
General Monitoring Cloth		●	●	●							●	
General Monitoring Foam		●	●	●							●	
General Monitoring Tape		●	●	●							●	
Radiotranslucent Cloth						●	●					
Radiotranslucent Foam						●	●					
Radiotranslucent Tape						●	●					
Stress/Diaphoretic Foam					●			●				
Radiotranslucent Stress/Diaphoretic Foam					●	●	●	●				
Wet Gel Stress/Diaphoretic Foam					●			●				●
Pediatric Foam									●			
Radiotranslucent Pediatric Foam						●	●		●			
Neonatal Electrodes										●		
Resting EKG Electrodes					●							

This list is not comprehensive and should be used as a recommended guideline only.

Point-of-Care Testing

Medline Point-of-Care testing products allow for quick monitoring, diagnoses and appropriate treatments for your patients without the inconvenience of waiting for results from the lab. These tests are fast, easy-to-use, accurate and manufactured under strict quality controls, so you get the reliability and accuracy your patients expect.



MPH1540

MPH1550



MPH3540

MPH3550

EvenCare® G2® Blood Glucose Monitoring System

Advanced technology makes diabetic care easier.

- » **Easy-to-use:** Voice assistance guides you through the testing process and provides results in either English or Spanish
- » **No coding required:** Minimizes errors and inaccurate results
- » **Fast test time:** Results in 6 seconds
- » **Small blood sample size:** 1 μ L
- » **Alternate site testing:** Option to test on the forearm, palm or fingertips
- » **Memory option:** Stores 300 test times and can provide 7-, 14- and 30-day averaging
- » **Alarm option:** Set reminders for up to 3 different times

Item No.	Description	Pkg.
MPH1540	Glucose Meter	1/ea
MPH1550	Test Strips	50/bx, 12 bx/cs
MPH1550Z	Test Strips	50/bx
MPH1560	Control Solution High/Low	1/bx, 6 bx/cs
MPH1560Z	Control Solution High/Low	1/bx
MPHBAAAZ	Replacement Battery, AAA	24/bx

EvenCare® G3 Professional Blood Glucose Monitoring System

The EvenCare G3 Professional Blood Glucose Monitoring System from Medline is designed to meet the unique needs of healthcare facilities and their patients. When used as directed, it's cleared for multi patient use.

- » **Quality control lock mode:** Ensures quality control testing is performed
- » **Test strip ejector:** Helps reduce cross contamination by eliminating contact with used strip
- » **No coding required:** Minimizes errors and inaccurate results
- » **Fast test time:** Results in 6 seconds
- » **Small blood sample size:** 0.8 μ L
- » **Alternate site testing:** Option to test on the forearm, palm or fingertips

Item No.	Description	Pkg.
MPH3540	Glucose Meter	1/ea
MPH3550	Test Strips	50/bx, 12 bx/cs
MPH3550Z	Test Strips	50/bx
MPH3550FL	Foil Wrapped Test Strips	50/bx, 12 bx/cs
MPH3550FLZ	Foil Wrapped Test Strips	50/bx
MPH3560	Control Solution High/Low	1/bx, 6 bx/cs
MPH3560Z	Control Solution High/Low	1/bx
MPHBKECR2032	Replacement Battery, 3V	6/bx

EvenCare Mini Meter

- » **Easy to use:** It provides the basic features to perform a glucose test
- » **Alternate site testing:** Can be used on the fingertip, palm, forearm, upper arm, calf, or thigh
- » **Before and after meal markers:** Mark your result according to when the test was performed; before eating or after eating.
- » **Discreet carrying pouch:** Fits in a purse or pocket

Item No.	Description	Pkg.
MPH5540	Glucose Meter	1/ea
MPH5550	Test Strips	50/bx, 20 bx/cs
MPH5550Z	Test Strips	50/bx
MPH5560	Control Solution Normal	1/bx, 50 bx/cs
MPH5560Z	Control Solution Normal	1/bx
MPHBKECR2032	Replacement Battery, 3V	6/bx



MPH5540

Coag-Sense™ PT/INR System

The Coag-Sense PT/INR (Prothrombin Time/International Normalized Ratio) system provides a more direct approach to INR testing as the only portable monitoring system that directly detects clotting endpoint, providing greater confidence in results. Features convenient testing with no coding required, rapid and reliable results and a wide hematocrit range of 15–60%.

Item No.	Description	Pkg.
ABO03P6001	Meter, Professional User's Manual, Professional Quick Reference Guide, Auto Single Use Lancets—starter bag of 12, 4 AA Alkaline batteries, AC power supply, carrying case	1/ea
ABO03PTESTKIT	Test Strip Kit—includes 54 transfer tubes (ABO03P5354Z)	50/bx
ABO03P6910Z	Control Strip Kit—includes 5 high and 5 low control strips and 1 vial of solution	10/bx
ABO03P5502Z	MiniPipette	2/pk
ABO03P5402Z	Pipette Tips	96/tr



ABO03P6001

Heel Incision Devices

- » High-quality steel blades make small, precise incisions for reduced bruising
- » Disposable, one-time use devices
- » Easy trigger activates blade quickly, safely and precisely
- » Color-coded for easy identification of devices


Item No.	Description	Pkg.
MPHMPREEMIE	Micro-Premie Heel Lancet, Light Blue	50/bx, 4 bx/cs
MPHPREEMIE	Premie Heel Lancet, White	50/bx, 4 bx/cs
MPHNEWBORN	Newborn Heel Lancet, Pink	50/bx, 4 bx/cs
MPHTODDLER	Toddler Heel Lancet, Purple	50/bx, 4 bx/cs



Monofilament

This monofilament can be used as a routine test to assess foot sensitivity for early signs of neuropathy.

Item No.	Description	Pkg.
MSC0002	Foot Monofilament, Diabetic Neuropathy Test	25/bx



General Purpose Lancets

Ultra sharp lancets have a fine-point tip for smooth and fast procedures. These lancets can be used with most universal lancing devices.

Item No.	Description	Pkg.
MPH100030	30 Gauge	100/bx, 25 bx/cs
MPH100030Z	30 Gauge	100/bx



Lancing Device

Fits all general purpose lancets.

Item No.	Description	Pkg.
MPHLANCINGPEN	Medline Adjustable Lancing Device	1/bx



Pen Needles

Medline pen needles provide easy medication delivery at every site.

- » Universal fit.
- » Thin wall technology allows for easy flow of medication during injection.

Item No.	Description	Pkg.
MPHPN2912	29 Gauge, 12 mm	100/bx, 10 bx/cs
MPHPN2912Z	29 Gauge, 12 mm	100/bx
MPHPN315	31 Gauge, 5 mm	100/bx, 10 bx/cs
MPHPN315Z	31 Gauge, 5 mm	100/bx
MPHPN318	31 Gauge, 8 mm	100/bx, 10 bx/cs
MPHPN318Z	31 Gauge, 8 mm	100/bx
MPHPN324	32 Gauge, 4 mm	100/bx, 10 bx/cs
MPHPN324Z	32 Gauge, 4 mm	100/bx



Safety Lancets

Medline safety lancets are simple and safe for both healthcare professionals and patients.

- » Needle retracts automatically and cannot be re-used protecting patients and healthcare professionals from needlestick injuries
- » Different sizes allow healthcare professionals to obtain the proper amount of blood for capillary testing
- » Quick and easy-to-use
- » Double spring action and controlled needle depth provides greater patient comfort

Item No.	Description	Pkg.
Push Button		
MPHSAFETY28	28 Gauge, 1.6 mm	200/bx, 18 bx/cs
MPHSAFETY28Z	28 Gauge, 1.6 mm	200/bx
MPHSAFETY285	28 Gauge, 1.6 mm	50/bx, 24 bx/cs
MPHSAFETY285Z	28 Gauge, 1.6 mm	50/bx
MPHSAFETY281	28 Gauge, 1.6 mm	100/bx, 20 bx/cs
MPHSAFETY281Z	28 Gauge, 1.6 mm	100/bx
MPHSAFETY21	21 Gauge, 1.8 mm	18 bx/cs
MPHSAFETY21Z	21 Gauge, 1.8 mm	200/bx

Pressure Activated

MPHST28	28 Gauge, 1.8 mm	10 bx/cs
MPHST28Z	28 Gauge, 1.8 mm	200/bx
MPHST21	21 Gauge, 2.2 mm	10 bx/cs
MPHST21Z	21 Gauge, 2.2 mm	200/bx

Side Button

MPHSTSIDE28	28 Gauge, 1.8 mm	200/bx, 10 bx/cs
MPHSTSIDE28Z	28 Gauge, 1.8 mm	200/bx
MPHSTSIDE21	21 Gauge, 2.4 mm	200/bx, 10 bx/cs
MPHSTSIDE21Z	21 Gauge, 2.4 mm	200/bx



Item No.	Description	Pkg.
Contact Activated		
MPHPRESS21	21 Gauge, 2.2 mm	100/bx, 20 bx/cs
MPHPRESS21Z	21 Gauge, 2.2 mm	100/bx
MPHPRESS23	23 Gauge, 2.2 mm	100/bx, 20 bx/cs
MPHPRESS23Z	23 Gauge, 2.2 mm	100/bx
MPHPRESS26	26 Gauge, 1.8 mm	100/bx, 20 bx/cs
MPHPRESS26Z	26 Gauge, 1.8 mm	100/bx
MPHPRESS28	28 Gauge, 1.8 mm	100/bx, 20 bx/cs
MPHPRESS28Z	28 Gauge, 1.8 mm	100/bx
MPHPRESS30	30 Gauge, 1.5 mm	100/bx, 20 bx/cs
MPHPRESS30Z	30 Gauge, 1.5 mm	100/bx
MPHPRESSBL	Lancet, Blade, 2.4 mm	100/bx, 20 bx/cs
MPHPRESSBLZ	Lancet, Blade, 2.4 mm	100/bx



AB004J6020



AB007P0501

i-STAT® System

Streamline the healthcare process from testing to diagnosis with the fully automated i-STAT system. A full diagnostic menu allows a single product to test for a variety of indicators to improve patient care. The small sample size of 2–3 drops of blood, the mobility and fast, lab quality results make the i-STAT a must have for every patient setting.

Item No.	Description	Pkg.
AB004J6020	i-STAT System Waived Kit	1/ea
AB009P3125	i-STAT Test Cartridge Chem 8+	25/bx
AB003P8425	i-STAT Test Cartridge Creatinine	25/bx
AB003P8325	i-STAT Test Cartridge Glucose	25/bx
AB003P8225	i-STAT Test Cartridge E3+	25/bx
AB003P8125	i-STAT Test Cartridge EC4+	25/bx
AB003P8025	i-STAT Test Cartridge 6+	25/bx

Piccolo Xpress® Chemistry Monitor

Reduce patient in-office time and improve office efficiency with the Piccolo Xpress Chemistry Monitor. Fast, three-step process provides lab quality results in a fraction of the time, partnered with an extensive testing menu, time with the patient can be maximized. The Piccolo Xpress system offers the most CLIA-Waived tests on a single platform.

Item No.	Description	Pkg.
AB007P0501	Piccolo Xpress Chemistry Monitor	1/ea
AB007P0208	Piccolo Reagent Disc Comprehensive Metabolic Panel	10/bx
AB007P0204	Piccolo Reagent Disc Basic Metabolic Panel	10/bx
AB007P0205	Piccolo Reagent Disc Lipid Panel	10/bx
AB007P0210	Piccolo Reagent Disc Liver Panel Plus	10/bx
AB007P0201	Piccolo Reagent Disc General Chem 6	10/bx
AB007P0207	Piccolo Reagent Disc Renal Panel	10/bx

**FOR MORE
CARTRIDGE
OPTIONS, PLEASE
CONTACT YOUR
MEDLINE REP.**

Influenza A & B Test

The Medline Influenza Test is a simple, rapid test that aids in the differential diagnosis of influenza A & B using nasal swab specimens. The easy-to-use format and rapid results allow for quick diagnosis and timely treatment.

Advantages

- » Results in 10 minutes
- » Dual reimbursement potential
- » Helps avoid unnecessary antibiotic use
- » Built-in control ensures accuracy
- » Sensitivity*
Flu A: 93.8% | Flu B: 77.4%
- » Specificity*
Flu A: 95.8% | Flu B: 98.0%
- » CLIA waived

Item No.	Description	Pkg.
MPH60025FLU	Medline Influenza A & B Test Strips	25/bx



Strep A Test

Medline Strep A Test Strips provide the rapid, qualitative detection of Strep A antigen using a throat swab. This allows for accelerated patient and resident care.

Advantages

- » Results in 5 minutes
- » Helps avoid unnecessary antibiotic use
- » Empowers timely treatment decisions
- » Built-in control ensures accuracy
- » Sensitivity*: 97%
- » Specificity*: 95%
- » CLIA waived

Item No.	Description	Pkg.
MPH30025SA	Medline Strep A Test Strips	25/bx



iFOB Tests

Simple and fast, this patient-friendly test is non-invasive, requires only one sample per test, and provides high sensitivity and specificity.

Advantages

- » Results in 5 minutes
- » Specific for human hemoglobin protein
- » No dietary or medication restrictions
- » Can be used for screening and diagnostic purposes
- » Built-in control ensures accuracy
- » 50ng hHB/mL analytical sensitivity*
- » 99.7% relative accuracy*
- » CLIA waived

Item No.	Description	Pkg.
MPH70025CTM	Test Cassettes, Collection Tubes and Mailers	25/bx
MPH70025CT	Test Cassettes and Collection Tubes	25/bx
MPH70025M	iFOB Mailers	25/bx
MPH70025T	iFOB Collection Tubes	25/bx



Pregnancy Tests

Medline provides a range of highly accurate, rapid tests that detect human chorionic gonadotropin (hCG) in urine or urine/serum for early detection of pregnancy.


Advantages

- » Results in 3 minutes (5 minutes for serum); positive may become visible sooner
- » Greater than 99% sensitivity*
- » Reduces staff hands-on time and increases office efficiency
- » Provides flexibility to meet your testing needs with a variety of options
- » Built-in controls ensure accuracy
- » hCG detection*:
 - Urine strip and cassette tests detect 25 mIU/ml
 - Combo+ test detects 20 mIU/ml for urine and 10 mIU/ml for serum
- » CLIA complexity:
 - Waived using urine
 - Moderate complexity using serum

Item No.	Description	Pkg.
MPH22025HCG	Medline hCG Pregnancy Combo+ Test Cassette	25/bx
MPH12025HCG	Medline hCG Pregnancy Test Cassette	25/bx
MPH11025HCG	Medline hCG Pregnancy Test Strip	25/bx



MPH22025HCG



MPH12025HCG



MPH11025HCG

*Refer to Instructions for Use

Medline 120 Urine Analyzer

The Medline 120 Urine Analyzer is specifically designed to improve accuracy and efficiency, allowing you and your staff to spend less time testing and more time with patients. Analyzer reading also helps standardize patient care.

Advantages

- » Automatic calibration for accurate results and easy operation
- » Reads Medline Urinalysis 40B, 50B, URO7 and 10SG parameter strips
- » Saves and recalls the last 2,000 results
- » Optional barcode reader
- » User and QC lockout

Item No.	Description	Pkg.
MPHUA120PROMO	Medline 120 Urine Analyzer Promo (includes analyzer and 5 bottles of 10SG strips)	1/kt



Medline Urinalysis Reagent Strips

The easy-to-read strips provide fast, qualitative results, aiding in the preliminary diagnosis and screening of many potential health issues.

Advantages


- » Multiple strip configurations to fit your testing needs
- » Compatible for visual and analyzer reading
- » Analytes appear in the same order as leading brands
- » CLIA waived

Item No.	Description	Pkg.
MPHUA10SG	Medline 10SG Urinalysis Strips	100/bt
MPHUAURO7	Medline URO7 Urinalysis Strips	100/bt
MPHUA50B	Medline 50B Urinalysis Strips	100/bt
MPHUA40B	Medline 40B Urinalysis Strips	100/bt
MPHUA2GP	Medline 2GP Urinalysis Strips (visual read only)	100/bt



Analytes

Item No.	Leukocytes	Nitrite	Urobilinogen	Protein	pH	Blood	Specific Gravity	Ketone	Bilirubin	Glucose
MPHUA10SG	•	•	•	•	•	•	•	•	•	•
MPHUAURO7	•	•	•	•	•	•		•		•
MPHUA50B	•	•		•		•				•
MPHUA40B	•	•		•						•
MPHUA2GP				•						•



Microscope Slides

- » Produced from premium glass sheets
- » Packaged in boxes of 72 each
- » Moisture absorbent sheet in each box to prevent sticking
- » Pre-cleaned and ready to use
- » Ground 90° edges
- » Dimensions: 25 mm x 75 mm x 1 mm

Item No.	Description	Pkg.
MLAB1301	Plain	1,440/cs, 20 bx/cs
MLAB1304	Frosted on One Side	1,440/cs, 20 bx/cs
MLAB1304W	White Frosted	1,440/cs, 20 bx/cs



Transfer Pipettes

- » Produced from premium glass sheets
- » Unbreakable
- » Packaged in a dispenser box, ideal for bench top use

Item No.	Description	Pkg.
MLAB138050	5 ml Non-Graduated, 2 ml Bulb Draw	500/bx, 10 bx/cs
MLAB138050Z	5 ml Non-Graduated, 2 ml Bulb Draw	500/bx
MLAB137040	5 ml Graduated, 2 ml Bulb Draw	500/bx, 10 bx/cs
MLAB137040Z	5 ml Graduated, 2 ml Bulb Draw	500/bx
MLAB135030	7 ml Graduated, 3.2 ml Bulb Draw	500/bx, 10 bx/cs
MLAB135030Z	7 ml Graduated, 3.2 ml Bulb Draw	500/bx



Cover Glass

- » Each flip top box contains a desiccant; vacuum-sealed in aluminum to prevent sticking
- » Packaged 1-oz. per box, 10 boxes per carton
- » #1 thickness (0.13 mm – 0.16 mm)

Item No.	Description	Pkg.
MLAB140410	24 x 50 mm	10/cs
MLAB141510H	24 x 50 mm	1/ea
MLAB140410H	22 x 22 mm	1/ea



Glass Test Tube

- » High quality borosilicate glass reduces pH changes and minimizes leaching concerns compared to soda lime glass
- » Neatly packaged in heavy duty boxes of 250 tubes each
- » Shrink wrapped to keep tubes particle-free
- » Compact boxes fit into most test tube dispensers
- » Well-rounded bottoms and smooth, fire-polished rims

Item No.	Description	Pkg.
MLAB1503	3 ml, 10 mm x 75 mm	1,000/cs, 4 bx/cs
MLAB1505	5 ml, 12 mm x 75 mm	1,000/cs, 4 bx/cs
MLAB1510	7 ml, 13 mm x 100 mm	1,000/cs, 4 bx/cs



Lab Marker

Item No.	Description	Pkg.
MLAB140020	Lab Marker	10/bx



DLS6285

15mL Centrifuge Tubes

- » Black printed graduations in 0.5 ml increments
- » White printed writing space
- » Polypropylene

Item No.	Description	Pkg.
DLS6285	Graduated 15 ml PP, Sterile	500/cs
DLS6286	Graduated 15 ml PP, Sterile	500/cs
DLS6284	Graduated 15 ml PP, Non-Sterile, Bulk	500/cs



DLS6288

50mL Centrifuge Tubes

- » Black printed graduations in 0.5 ml increments
- » White printed writing space
- » Polypropylene

Item No.	Description	Pkg.
DLS6288	Graduated 50 ml PP, Sterile	500/cs
DLS6289	Graduated 50 ml PP, Sterile	500/cs
DLS6287	50 ml, Non-Sterile, Bulk	500/cs



Biohazard Specimen Bags

Clear specimen bags with zip closure, dual pouch and biohazard symbol.

Item No.	Description	Pkg.
DYND30260	6" x 6" Bag w/ 6" x 6" Document Pouch	1,000/cs
DYND30261	6" x 9" Bag w/ 6" x 9" Document Pouch	1,000/cs
DYND30271	8" x 8" Bag w/ 8" x 8" Document Pouch	1,000/cs
DYND30281	12" x 15" Bag w/ 12" x 15" Document Pouch	500/cs



Pathology Containers with Lids

- » Made of high-grade, durable polypropylene
- » May be frozen without cracking
- » Secure, snap-lock caps minimize leaking

Item No.	Description	Pkg.
DYND34255	8-oz.	250/cs
DYND34260	16-oz.	100/cs
DYND34265	32-oz.	100/cs
DYND34270	64-oz.	50/cs
DYND34275	165-oz.	25/cs



DYND30382

DYND30330

DYND30370

PMR30363

DYND30362

DYND30100

Click-N-Close Pneumatic Tube-Safe Specimen Containers

The continuous-thread lid on these 95kPA specimen containers twist close with an audible click to ensure a tight and secure fit.

Item No.	Description	Pkg.
DYND30382	3-oz. (90 ml), Sterile, Tamper-evident Seal	400/cs
DYND30387	3-oz. (90 ml), Individually Packaged in Peel Pouch, Appropriate for O.R. Use	100/cs
DYND30385	4-oz. (120 ml), Sterile, Tamper-evident Seal	300/cs
DYND30389	4-oz. (120 ml), Individually Packaged in Peel Pouch, Appropriate for O.R. Use	100/cs

Dual Click Pneumatic Tube-Safe Specimen Containers

A distinct double click lets you know that these 95kPa specimen containers are securely sealed and ready for transport.

Item No.	Description	Pkg.
PMR30363	3-oz. w/Tamper-evident Seal, Sterile	400/cs
PMR30364	4-oz. w/Tamper-evident Seal, Sterile	300/cs

General Use Specimen Containers

These economical containers are easy to use and meet a wide variety of specimen collection needs.

Item No.	Description	Pkg.
DYND30330	4-oz. (120 ml), Individually Packaged in Poly Bag, Sterile, ID Label	100/cs
DYND30335	4-oz. (120 ml), Container and Lid Packaged Separately in Sleeves of 25, Non-Sterile	500/cs
DYND30350	4.5-oz. (130 ml), Individually Packaged in Poly Bag, Sterile, ID Label	100/cs

Pneumatic Tube-Safe Specimen Containers

Ideal for pneumatic tube systems, our 95kPA-compliant specimen containers feature a continuous threaded screw-on lid for a tight, leak-resistant seal.

Item No.	Description	Pkg.
DYND30362	3-oz. Sterile	400/cs
DYND30367	3-oz., O.R. Sterile, Packaged Individually in Peel Pouch	100/cs
DYND30364	4-oz. Sterile	300/cs
DYND30369	4-oz., O.R. Sterile, Packaged Individually in Peel Pouch	100/cs

Wide-Mouth Click-N-Close Pneumatic Tube-Safe Specimen Container

Our wide-mouth 95kPa specimen containers seal tightly with a noticeable click, confirming that the lid is tightly sealed against leaks.

Item No.	Description	Pkg.
DYND30370	100 ml, Sterile Fluid Pathway with Tamper-evident Patient ID Band	400/cs

Non-Sterile Urinalysis Containers

A convenient pouring spout and a broad base helps minimize accidental spilling and tipping.

Item No.	Description	Pkg.
DYND30100	6-oz., w/o Lid	400/cs
DYND30102	6-oz., w/ Pre-Printed Lids	100/cs
DYND30103	6-oz., w/ Pre-Printed Lids	300/cs



Medline Industries, Inc. One Medline Place, Mundelein, IL 60060

FOLLOW US [f](#) [in](#) [t](#) [BLOG](#)

Medline United States
1-800-MEDLINE (633-5463)
medline.com | info@medline.com

Medline Canada
1-800-396-6996
medline.ca | canada@medline.com

Medline México
01-800-831-0898
medlinemexico.com | mexico@medline.com

We reserve the right to correct any errors that may occur within this brochure. Some products may not be available for sale in Mexico or Canada. ©2012, 2013, 2014, 2016 Medline Industries, Inc. All rights reserved. Agilent is a registered trademark of Agilent Technologies, Inc. BCI is a registered trademark of BCI, Inc. CAS and Casmed are registered trademarks of Cas Medical Systems, Inc. Coag-Sense is a trademark of Coagusense, Inc. Colin is a registered trademark of Omron Healthcare Co. Criticare is a registered trademark of Criticare Systems, Inc. Critikon, Dinamap and Marquette are registered trademarks of GE Medical Systems Information Technologies, Inc. Datex-Ohmeda is a registered trademark of GE Healthcare Finland Oy Corporation. Defibtech, Lifeline AED and Lifeline VIEW are registered trademarks of Defibtech, LLC. Dräger is a registered trademark of Draegerwerk AG & Co. E-scope is a registered trademark of Excellus Technologies, Inc. Fukuda Denshi is a registered trademark of Fukuda Denshi Co. GE is a registered trademark of General Electric Company. i-STAT is a registered trademark of the Abbott Group of Companies. Nihon Kohden is a registered trademark of Nihon Kohden Corporation. Philips is a registered trademark of Koninklijke Philips Electronics N.V. Limited Liability Company. Piccolo Xpress is a registered trademark of Abaxis, Inc. R.A.T.E. is a trademark of BestMed. Siemens is a registered trademark of Siemens Aktiengesellschaft. Spacelabs is a registered trademark of Spacelabs Healthcare. TemporalScanner is a trademark and Exergen is a registered trademark of Exergen Corporation. Tyvek is a registered trademark of E. I. Du Pont De Nemours And Company. Welch Allyn is a registered trademark of Welch Allyn, Inc. Complimates, Evencare, G2, MedGel, Medline and Synergy Classic are registered trademarks of Medline Industries, Inc. MKT1664873 / LIT007A / 3M / K&M / 3