

# Suction Catheters, Tracheostomy Care and Cleaning Trays

Essential products for effective airway clearance.



**Clear the  
way to  
better care**

# Meeting all your airway management needs

Patients with artificial airways and tracheostomies are at high risk for airway obstruction, impaired ventilation and infection. Skilled bedside nursing care can help prevent these complications.<sup>1</sup>

We can support clinicians with a comprehensive portfolio of suction catheters in a wide range of designs, sizes and kit configurations. We also offer a variety of tracheostomy care and cleaning trays to satisfy clinician preference and simplify supply management.



## Customize your trays for maximum efficiency

Save time, reduce waste and stretch your budget dollars.



## Leverage stock kits and trays

Streamline supply management and ordering. Stock kits and trays provide single-item convenience and best practice standardization.



## Reduce costs

As the #1 supplier of custom trays, we manage the entire supply chain—from sourcing and manufacturing to warehousing and distribution. This means stability and lower prices for our customers.

**What's  
inside**

<b>Open Suction Catheters .....</b>	<b>5</b>
<b>Sleeved Catheters .....</b>	<b>5</b>
<b>Kits and Trays .....</b>	<b>6</b>
<b>Suction Catheter Kits with Sterile Solutions .....</b>	<b>7</b>
<b>Addipak® Sterile Solutions.....</b>	<b>8</b>
<b>Irrigation Solutions .....</b>	<b>8</b>
<b>Vac-Assist Portable Aspirator.....</b>	<b>9</b>
<b>Vacu-Line Suction Aspirator .....</b>	<b>9</b>
<b>Mucus/Specimen Traps .....</b>	<b>10</b>
<b>Suction Aspirator .....</b>	<b>10</b>
<b>Premium Tracheostomy Tube Holders.....</b>	<b>11</b>
<b>Standard Tracheostomy Tube Holders .....</b>	<b>11</b>
<b>Tracheostomy Care and Cleaning Trays.....</b>	<b>12</b>
<b>Tracheostomy Care Components Checklist.....</b>	<b>13</b>

# Contro-Vac valve helps reduce risk of infection

## Better design drives better performance

All of our tracheal suction catheters feature our exclusive Contro-Vac valve to help reduce the risk of infection, airway trauma and patient discomfort.

With its 5-in-1 connector and stabilizing bar, the Contro-Vac valve connects easily and securely to any size tubing. Forward placement of the valve reduces risk of cross-contamination from fluids.

## Proven to deliver enhanced protection

Independent testing compared Medline's Contro-Vac valve to the leading catheter's chimney valve. Contro-Vac was more effective at preventing the release of contaminants during suctioning.<sup>2</sup>



**Medline Contro-Vac valve**



**Chimney valve from leading competitor**



# Open Suction Catheters


- Available in a complete range of adult and pediatric sizes
- All Medline suction catheters feature easy-to-read depth markings
- DeLee tip is flat, soft and smooth to minimize airway trauma for newborns, infants and children
- Angled, soft and smooth whistle tip maximizes suctioning capability for adult patients, while minimizing trauma potential

## Single Sterile-Straight Pack

Item No.	Description	Pkg.
DYND41906	Suction Catheter, Single Sterile: 6 Fr	100/cs
DYND41908	Suction Catheter, Single Sterile: 8 Fr	100/cs
DYND41900	Suction Catheter, Single Sterile: 10 Fr	100/cs
DYND41901	Suction Catheter, Single Sterile: 12 Fr	100/cs
DYND41902	Suction Catheter, Single Sterile: 14 Fr	100/cs
DYND41903	Suction Catheter, Single Sterile: 16 Fr	100/cs
DYND41904	Suction Catheter, Single Sterile: 18 Fr	100/cs
DYND40990	Red Poly-Cath Open Suction Catheter Straight Pack: 10 Fr	100/cs
DYND40992	Red Poly-Cath Open Suction Catheter Straight Pack: 14 Fr	100/cs



DYND41902



DYND41900



DeLee Tip



Whistle Tip

Medline was the first in the industry to offer depth mark catheters.

# Sleeved Catheters

Our D407 series open-suction catheters are contained in a sleeve to help avoid cross-contamination.

Item No.	Description	Pkg.
DYND40708F	Sleeved Suction Catheter: 8 Fr	50/cs
DYND40700F	Sleeved Suction Catheter: 10 Fr	50/cs
DYND40701F	Sleeved Suction Catheter: 12 Fr	50/cs
DYND40702F	Sleeved Suction Catheter: 14 Fr	50/cs



DYND40701F

## Kits and Trays

- Includes all the necessities for airway clearance
- Irrigation solutions can be poured into the tray or the pop-up Cup
- All kits and trays contain Medline Aloetouch® 3G Gloves
- Catheters have easy-to-read depth markings
- DeLee tip is flat, soft and smooth to minimize airway trauma for newborns, infants and children
- Angled, soft and smooth whistle tip maximizes suctioning capability for adult patients, while minimizing trauma potential

Medline is the world's largest kit manufacturer—making more than 600,000 surgical and minor procedure trays every day.

### Suction Catheter Trays

Item No.	Description	Pkg.
DYND40986	Suction Catheter Mini Tray with Gloves: 6 Fr	100/cs
DYND40988	Suction Catheter Mini Tray with Gloves: 8 Fr	100/cs
DYND40980	Suction Catheter Mini Tray with Gloves: 10 Fr	100/cs
DYND40981	Suction Catheter Mini Tray with Gloves: 12 Fr	100/cs
DYND40982	Suction Catheter Mini Tray with Gloves: 14 Fr	100/cs
DYND40983	Suction Catheter Mini Tray with Gloves: 16 Fr	100/cs
DYND40984	Suction Catheter Mini Tray with Gloves: 18 Fr	100/cs



### Suction Catheter Kits

Item No.	Description	Pkg.
DYND40976	Suction Catheter Pack with Gloves and Pop-up Cup, 6 Fr	50/cs
DYND40978	Suction Catheter Pack with Gloves and Pop-up Cup, 8 Fr	50/cs
DYND40970	Suction Catheter Pack with Gloves and Pop-up Cup, 10 Fr	50/cs
DYND40971	Suction Catheter Pack with Gloves and Pop-up Cup, 12 Fr	50/cs
DYND40972	Suction Catheter Pack with Gloves and Pop-up Cup, 14 Fr	50/cs
DYND40973	Suction Catheter Pack with Gloves and Pop-up Cup, 16 Fr	50/cs
DYND40974	Suction Catheter Pack with Gloves and Pop-up Cup, 18 Fr	50/cs



### Suction Catheter Trays, Poly-cath

Item No.	Description	Pkg.
DYND48980A	Suction Catheter Mini Tray with Gloves, Poly-Cath, 10 Fr	100/cs
DYND48982	Suction Catheter Mini Tray with Gloves, Poly-Cath: 14 Fr	100/cs



# Suction Catheter Kits with Sterile Solutions

- Kits include everything you need for airway clearance
- All kits and trays contain Medline Aloetouch 3G Gloves
- Catheters have easy-to-read depth markings
- DeLee tip is flat, soft and smooth to minimize airway trauma for newborns, infants and children
- Angled, soft and smooth whistle tip maximizes suctioning capability for adult patients, while minimizing trauma potential

## Water

Item No.	Description	Pkg.
DYND41476	Suction Catheter Tray with Sterile Water, 100 mL, 6 Fr	36/cs
DYND41478	Suction Catheter Tray with Sterile Water, 100 mL, 8 Fr	36/cs
DYND41470	Suction Catheter Tray with Sterile Water, 100 mL, 10 Fr	36/cs
DYND41471	Suction Catheter Tray with Sterile Water, 100 mL, 12 Fr	36/cs
DYND41472	Suction Catheter Tray with Sterile Water, 100 mL, 14 Fr	36/cs
DYND41474	Suction Catheter Tray with Sterile Water, 100 mL, 18 Fr	36/cs



## Saline

Item No.	Description	Pkg.
DYND41478SA	Suction Catheter Tray with Sterile Saline, 100 mL, 8 Fr	36/cs
DYND41470SA	Suction Catheter Tray with Sterile Saline, 100 mL, 10 Fr	36/cs
DYND41442	Suction Catheter Tray with Sterile Saline, 100 mL, 14 Fr	36/cs
DYND40442	Suction Catheter Tray with Saline and Drape, 14 Fr	36/cs



# Addipak® Sterile Solutions

- Tamper-proof, color-coded vials for easy solution identification
- Inverted mL graduations and precise opening for accurate and controlled dispensing solution
- 15 mL vials packaged 48 to a box, 3 boxes to a case
- 3 mL and 5 mL vials packaged 100 to a box, 10 boxes to case
- Not for use in parenteral applications

Item No.	Description	Pkg.
HUD20034	0.45% Sterile Saline (Green), 3 mL	1,000/cs
HUDRHU39U	0.9% Sterile Saline (Pink), 3 mL	1,000/cs
HUD20031	Sterile Water (Blue), 3 mL	1,000/cs
HUDRHU59U	0.9% Sterile Saline (Pink), 5 mL	1,000/cs
HUD20051	Sterile Water (Blue), 5 mL	1,000/cs
HUDRHU69U	0.9% Sterile Saline (Pink), 15 mL	144/cs



Sterile Water



0.9% Sterile Saline



0.45% Sterile Saline

# Irrigation Solutions

- Ready-to-use containers reduce the contamination risks and potential waste associated with bulk containers of sterile water and saline
- Ideal for device irrigation, suction, dressing moistening and wound debridement
- Cup size provides the proper amount for suctioning

Item No.	Description	Pkg.
RDI30295	Sterile Water, 100 mL	48/cs
RDI30296	Sterile Saline, 0.9%, 100 mL	48/cs



RDI30296



## Vac-Assist Portable Suction Aspirator

- Durable and portable suction aspirator designed specifically for home care and nursing homes
- Suction Range: up to 74.6 ± 3% kPa
- Pumping Rate: > 40 L/min
- Includes 800 mL/cc canister
- Operates on AC 120 V/60 Hz power and does not have a battery; includes medical grade cord



Item No.	Description	Pkg.
HCS7000	Vac-Assist Suction Aspirator with 800 mL Canister	1 ea

## Vacu-Line Suction Aspirator

- Updated, lightweight, compact aspirator ideal for home care
- Suction Range: 21.15 ± 4.86 to 82.84 ± 1.84 kPa
- Pumping Rate: = 14 L/min
- Canister Size: 850 mL
- Operates on AC 120 V/60 Hz power and does not have a battery; includes non-medical grade cord



Item No.	Description	Pkg.
HCS7001	Vacu-Line Suction Aspirator with 850 mL Canister	1 ea

## Suction Aspirator Accessories

- Easy-to-use, ultra-quiet aspirators come with a full assortment of accessories for diverse applications
- Weight: 12.1 lb./5.49 kg
- Maximum suction: -21 inHg/520 mmHg
- Rugged construction is ideal for bedside carts

Item No.	Description	Pkg.
HCS7850	Canister with Float Lid, 850 cc	12/cs
HCS7851	Canister with Float Lid and Tubing, 850 cc	12/cs
HCS7852	Canister with Filter Lid and Tubing, 850 cc	50/cs
HCS7853	Canister with Filter Lid, 850 cc	100/cs
HCS7120	Canister with Float Lid, 1,200 cc	12/cs
HCS7121	Canister with Float Lid and Tubing, 1,200 cc	12/cs



HCS7850


HCS7120



Tubing included with HCS7851, HCS7852 and HCS7121.

# Mucus/Specimen Traps

- Leak-proof transport caps for more secure specimen collection
- Sterile plastic vials, graduated in 1cc increments
- Includes threaded cap, tubing not made from natural rubber latex and patient ID label




DYND44108

## Sterile DeLee Mucus Trap

Item No.	Description	Pkg.
DYND44108	Trap with Catheter and Cap, Capacity 20 cc, 8 Fr	50/cs
DYND44110	Trap with Catheter and Cap, Capacity 20 cc, 10 Fr	50/cs

## Sterile Mucus Specimen Trap

Item No.	Description	Pkg.
DYND44140	Trap with Transport Cap, Capacity: 40 cc	50/cs
DYND44180	Trap with Transport Cap, Capacity: 80 cc	50/cs



DYND44180

# Premium Tracheostomy Tube Holders


- Hypoallergenic, no-stretch cotton straps help keep tracheostomy tubes securely in place for uninterrupted airflow
- Adjustable hook-and-loop closures
- Soft cotton helps ensure patient comfort while reducing skin irritation

## One-piece

Item No.	Description	Pkg.
HCS597LA	Large Adult, Fits 13–19" Neck	10/bx
HCS597M	Small Adult, Fits 9–14" Neck	10/bx
HCS597S	Neonatal/Pediatric, Fits 7–9" Neck	10/bx

## Two-piece

Item No.	Description	Pkg.
HCS540	Adult, Fits up to 20" Neck	10/bx
HCS541	Pediatric, Fits up to 18" Neck	10/bx
HCS542	Neonatal/Infant, Fits up to 9" Neck	10/bx



HCS540





HCS597S

# Standard Tracheostomy Tube Holders

- Soft, comfortable straps have adjustable hook-and-loop closures to keep tracheostomy tubes securely in place, ensuring uninterrupted airflow
- Available in two sizes to fit infants, young children and adults

Item No.	Description	Pkg.
ZPP301P	Trach Tube Holder, 1-Piece Pediatric	20/cs
ZPP301	Trach Tube Holder, 1-Piece Adult	20/cs
ZPP501P	Trach Tube Holder, 2-Piece Pediatric	20/cs
ZPP501	Trach Tube Holder, 2-Piece Adult	20/cs
ZPP401	Anti-Disconnect Endotracheal Tube Holder	20/cs



# Tracheostomy Care and Cleaning Trays

- Medline offers a wide variety of adult and pediatric tracheostomy care and cleaning trays in standard and customized configurations to meet diverse clinical needs.
- Available in economy, basic and full-compartment kits
- Ask about customization

	DYND40582	DYND40622	DYND40610	DYND40620	DYND40589	DYND40615	DYND40580	DYND40650	DYND40511	DYND40594	CC3T3030A	CC3T4691A	CC4682A	CC4681A
Tray Compartments	3	1	1	1	3	3	3	3	3	3	1	3	3	3
AloeTouch Gloves	2		2	2	2	2	2	2	2	2		2		2
Drape	1		1	1	1	1	1		1	1		1	1	1
Pipe Cleaners		2	2							3	3	3	3	3
Non-Woven Gauze 4" x 4"	4	2	2	2	2	2	2	2	1	3	4	1	3	3
Cotton Tip Applicator	2	2	5	5	2	5	2	2	4	2	2	2	3	3
Pliable Cleaning Brush	1	1	1		1	1	1		1	1	1	1	1	1
Twill Tape, 34"	1	1	1	1	1	1	1	1	1		1	1	1	1
6-Ply Trach Dressing	1		1	1	1	1	1	1	1	1	1	1	1	1
Hydrogen Peroxide					2	2		1		1				
Saline 100 mL					1				1	1				
Water 100 mL								1						
Suction Catheter 14 fr							1			1				
Forceps													1	1
Pkg.	20/cs	30/cs	20/cs	20/cs	20/cs	20/cs	20/cs	20/cs	20/cs	20/cs	30/cs	20/cs	20/cs	20/cs



# Tracheostomy Care Components Checklist

To create your own custom kit, choose among the components listed below. For additional help in designing your kit, contact your Medline Representative.

Account Name \_\_\_\_\_

Account Number \_\_\_\_\_



## Gloves



Description	Qty
Aloetouch® Vinyl 3G P/F Medium	
Vinyl P/F Medium	
Vinyl P/F Large	
Nitrile P/F Large	
Nitrile P/F Medium	

## Utensils

Description	Qty
Trach Brush	
6" Cotton Tip, Plastic	
4" Cotton Tip, Plastic	
Pipe Cleaners	

## Solutions

Description	Qty
100 mL Saline Pour Bottle	
100 mL Water Pour Bottle	
Hydrogen Peroxide, 1 oz	



## Catheters

Description	Qty
6 Fr Catheter	
8 Fr Catheter	
10 Fr Catheter	
12 Fr Catheter	
14 Fr Catheter	
16 Fr Catheter	

## Accessories

Description	Qty
Plastic Forcep	
Drape	
Towel, 13" x 19" 2PP, White	
Twill Tape, 1/4" x 24"	
Twill Tape, 1/2" x 36"	
CHLORASCRUB Swabstick, 1.6 mL	

## Drain Sponges

Description	Qty
Premium 4" x 4" Fenestrated Avant Gauze® Drain Sponge	
Premium 2" x 2" Fenestrated Avant Gauze Drain Sponge	
100% Cotton 4" x 4" Drain Sponge	



## Deluxe Avant Non-Woven Sponges

Description	Qty
2" x 2", 3-ply	
2" x 2", 4-ply	
3" x 3", 4-ply	
4" x 4", 3-ply	
4" x 4", 4-ply	

## Avant Non-Woven Sponges

Description	Qty
2" x 2", 3-ply	
2" x 2", 4-ply	
3" x 3", 4-ply	
4" x 4", 3-ply	
4" x 4", 4-ply	
4" x 4", 6-ply	

## Combos

Description	Qty
Set 1 Twill, 1 Drain Sponge, 2 4" x 4" Gauze	
Set 1 Twill, 1 Drain Sponge, 4 4" x 4" Gauze	

Additional components:







# Why customize your trays?



## Provide choice

Satisfy clinician preference and budget goals



## Create efficiencies

Eliminate time wasted picking additional items



## Increase value

Pay for only the products you know you'll use



## Reduce waste

Avoid discarding unused components

**Start today.** For more information and a free trial, contact your Medline Representative, visit [medline.com](http://medline.com) or call 1-800-MEDLINE.

**References:** 1. Nance-Floyd B. Tracheostomy care: An evidence-based guide to suctioning and dressing changes. American Nurse Today. 2011; 6(7). Available at: <http://www.americannursetoday.com/tracheostomy-care-an-evidence-based-guide-to-suctioning-and-dressing-changes/> Accessed April 29, 2015.  
2. Geneva Lab Study (1995). Data on file.



**We make  
healthcare  
run better™**

Medline Industries, LP | 1-800-MEDLINE (633-5463) | [www.medline.com](http://www.medline.com)

© 2024 Medline Industries, LP. All rights reserved. Medline is a registered trademark of Medline Industries, LP. All product names, regardless of whether or not they appear with a trademark symbol, are the property of their respective companies unless otherwise noted.  
MKT19W4978086 / LIT121R / 1.5M / SG / 87