

Advanced Wound Care Products

A wide assortment of options for HME providers



**Expand
your
offerings**

Advanced wound care products for HME providers

Our wide assortment of wound care products can help you expand your revenue stream and grow your business as the market rapidly shifts toward more patients receiving care at home.

As an innovative leader in advanced wound care products, we can help you build standardization of care with our unparalleled resources, breadth of products and expertise. This added value can help you attract new customers, build loyalty and bolster your cash business.

The U.S. home healthcare market is forecasted to grow **5% annually from 2020-2024.**¹

Five-layer silicone foam dressing

Secure yet gentle adherence to skin, with a foam core that ensures reliable absorption in wounds.

Item No.	Description	Pkg.
EQX2333	Five-layer silicone foam, 3" x 3" (7.5 x 7.5 cm)	10/bx, 100/cs
EQX2344	Five-layer silicone foam, 4" x 4" (10 x 10 cm)	10/bx, 100/cs
EQX2366	Five-layer silicone foam, 6" x 6" (15.2 x 15.2 cm)	10/bx, 100/cs
EQX2377	Five-layer silicone foam, Sacrum, 7" x 7" (17.8 x 17.8 cm)	5/bx, 40/cs
EQX2399	Five-layer silicone foam, Sacrum, 9" x 9" (22.8 x 22.8 cm)	5/bx, 25/cs



Silver alginate wound dressing

Absorbent alginate primary dressing for partial- and full-thickness wounds.

Item No.	Description	Pkg.
EQX9012	Silver alginate, 12" ribbon (30 cm)	5/bx, 50/cs
EQX9022	Silver alginate, 2" x 2" (5 x 5 cm)	10/bx, 100/cs
EQX9045	Silver alginate, 4" x 5" (10 x 12.7 cm)	10/bx, 100/cs
EQX9048	Silver alginate, 4" x 8" (10 x 20.3 cm)	5/bx, 50/cs
EQX9066	Silver alginate, 6" x 6" (15 x 15 cm)	5/bx, 50/cs



Alginate wound dressing

Absorbent alginate primary dressing for moderate to heavily draining partial- and full-thickness wounds.

Item No.	Description	Pkg.
EQX7012	Alginate, 12" ribbon (30 cm)	10/bx, 100/cs
EQX7022	Alginate, 2" x 2" (5 x 5 cm)	10/bx, 100/cs
EQX7044	Alginate, 4" x 4" (10 x 10 cm)	10/bx, 100/cs



Multi-layer compression system

Delivers therapeutic compressions for the treatment of chronic venous insufficiency and associated edema.

Item No.	Description	Pkg.
EQX4300	Compression, 3-layer	1/bx, 12/cs
EQX4400	Compression, 4-layer	1/bx, 12/cs



Bordered gauze dressing

Absorbent island dressing with a soaker pad and conformable, non-woven fabric with an adhesive backing for use as a primary or secondary dressing.

Item No.	Description	Pkg.
EQX3223	Bordered gauze, 2" x 3.5" (5 x 9 cm)	50/bx, 500/cs
EQX3244	Bordered gauze, 4" x 4" (10 x 10 cm)	25/bx, 500/cs
EQX3266	Bordered gauze, 6" x 6" (15 x 15 cm)	25/bx, 500/cs



Wound cleanser

Formulated with a gentle surfactant to allow fast, thorough and gentle cleansing of dermal wounds. Spray delivery loosens wound debris and is noncytotoxic, so it won't hurt healthy cells.

Item No.	Description	Pkg.
EQX6008	Wound cleanser, 8oz. bottle with trigger sprayer	6/cs



Hydrogel dressing

Gel-based wound care dressing donates moisture to help create an optimal wound healing environment.

Item No.	Description	Pkg.
EQX6102	Hydrogel, 1-oz. tube (28g)	30/cs



Learn more. Contact your Medline representative, call 1-800-MEDLINE or visit [medline.com](https://www.medline.com).



**We make
healthcare
run better.™**

Medline Industries, Inc.
Three Lakes Drive, Northfield, IL 60093
Medline United States | 1-800-MEDLINE (633-5463)
[medline.com](https://www.medline.com) | info@medline.com

Medline Canada
1-800-268-2848 | [medline.ca](https://www.medline.ca) | canada@medline.com
Medline México
01-800-831-0898 | [medlinemexico.com](https://www.medlinemexico.com) | mexico@medline.com

Follow us    

Reference: 1. PR Newswire. US Home Healthcare Market (Services, Equipment & Software): Insights & Forecast with Potential Impact of COVID-19 (2020-2024). May 4, 2020. Available at: <https://www.prnewswire.com/news-releases/us-home-healthcare-market-services-equipment--software-insights--forecast-with-potential-impact-of-covid-19-2020-2024-301052261.html>. Accessed November 25, 2020.

Some products may not be available for sale in Mexico or Canada. We reserve the right to correct any errors that may occur within this brochure.
© 2020 Medline Industries, Inc. All rights reserved. Medline is a registered trademark of Medline Industries, Inc. MKT19W1064131 / e20847 / 3